

Sediment trap, Ø110/104 mm 28.300 series

Manual



Research Equipment Limnology • Oceanography • Hydrobiology

	Manual for Ø110 mm sediment traps	Series 28.XXX
	Caution Do not use any kind of alcohol for cleaning the tubes, a mild soap is recommended. KC Denmark A/S is not, and cannot be held, responsible for any damage(s) made to equipment or to operators who ignore safety precautions or because of misuse or wrong operation.	
	Preparing the sediment trap	
1	The sediment trap can be configured in various ways, please see below for the requirements and the corresponding order numbers:	
2	 1 main rack 28.303 1 arm with pivot 28.304 1 sediment trap 28.301 1 steering fin, large 28.305 8 unbraco screws, M6 x 12, 28.109 	
3	 1 main rack 28.303 2 arms with pivot 28.304 2 sediment traps 28.301 8 unbraco screws, M6 x 12, 28.109 	

4	 1 main rack 28.303 2 arms with pivot 28.304 2 sediment traps 28.301 1 steering fin, small 28.305 A 1 balance arm 28.307 16 unbraco screws, M6 x 12, 28.109 	
5	 1 main rack 28.303 3 arms with pivot 28.304 3 sediment traps 28.301 1 steering fin 28.301 1 steering fin 28.305 16 unbraco screws, M6 x 12, 28.109 	
6	 1 main rack 28.303 4 arms with pivot 28.304 4 sediment traps, 28.301 16 unbraco screws, M6 x 12, 28.109 	

7	Mount all accessories to the main rack; it is highly recommended to locate the zinc anode at top. The zinc anode will protect the traps from corrosion.	
8	Attach the requested accessories by means of the 6 x 12 mm bolts; tighten the bolts "B, C, D and E" firmly. It is highly recommended that you add some antiseize compound like Never-Seez, copper paste or a similar product to the bolts thread. It will prevent the bolts from being stuck and you can easily remove these when changing the configuration. "A" shows the zinc anode. The zinc anode will corrode during deployment, and you must check that the zinc anode is intact and has a reasonable size. Furthermore, you may consider if the material left will be adequate for the time of deployment.	A B C D E
9	The pivot arm and the steering fin are attached on the main rack. You must choose the correct steering fin (2 sizes available) and the corresponding counterpart to keep the balance of the trap. The photo shows a smaller 80 mm version but basic mounting is equal.	

10	Before sliding in the tubes, unscrew the American wing nut, "F", so the bracket "G" becomes loose.	
11	While inserting the tube, hold the bracket "G" out of the tube, then, push back "G" so it resides on the collar "H" of the tube. Fasten the American wing nut "F" to secure the sample tube.	F H
12	Slide all tubes into the nylon holders and check if the pivots are operating freely. Check that all bolts and the American wing nuts "F" are fastened properly.	

13	Depending of your project you may apply some preservatives to each tube. Use a syringe with a silicone tube or similar. Before deployment, you may slowly fill each tube with water so the preservatives remain undisturbed at the bottom of the tube.	
	Various ways of deployment	
14	You can deploy a single sediment trap or you may connect more traps in cascade to have samples from various depths. For long terms sampling we will recommend one of the below-mentioned options:	
15	For a well-defined depth, you can attach an appropriate weight at the end of the line, and the buoy will keep the deployment line in vertical position.	Anchored mooring Buoy Pycnocline Weight

		Drifting mooring Buoy
16	For larger depth, we will recommend this option.	Pycnocline
		U Weight
	Emptying the traps	
17	The sediment sample can now be washed out with fresh water or distilled liquids. Caution Do not use any kind of alcohol for cleaning the tubes; a mild soap is recommended.	
	Maintenance	
18	Caution Do not use any kind of alcohol for cleaning the tubes, a mild soap is recommended. Before storing the sediment trap, it is	
Poved	important to clean all parts with fresh water, especially if the trap has been used in salt water. Remove the sample tube from its holder. Do not expose the sample tube to direct sunlight for a long period. March 21, 2019 – Ikj	

Rev.: March 21, 2019 - lkj



Research Equipment
Limnology • Oceanography • Hydrobiology

E-mail: kc@kc-denmark.dk website: http://www.kc-denmark.dk/
Holmbladsvej 17-19, DK 8600 Silkeborg. Denmark. Tel. +45 86 82 83 47. Fax +45 86 82 49 50
Bank: Sydbank. SWIFT: SYBKDK22 IBAN DK5070460000104832
VAT no. DK 29 61 96 62