

# Magnetic Stirrer PCE-MSR 50



## Mini Magnetic Stirrer PCE-MSR 50

**Mini magnetic stirrer to 1000 ml / 88 mm platform / Mains operation /  
Adjustable stirring speed**

The mini magnetic stirrer PCE-MSR 50 is suitable for applications in small laboratories. The mini magnetic stirrer can safely mix stirring quantities of 100 ... 1000 ml. Vessels with a diameter of up to 88 mm can be placed on the stirring surface of the mini magnetic stirrer PCE-MSR 50.

So the magnetic stirrer can be used for all common beakers. The magnetic stirrer is supplied with an external 5 V mains adapter. Due to its small size and weight, the mini magnetic stirrer can also be used in mobile field surveys.

- ▶ Magnetic stirrer up to 1000 ml
- ▶ Small dimensions
- ▶ Mains operation
- ▶ Adjustable speed to 1600 rpm
- ▶ Stirring surface 88 mm
- ▶ For field and laboratory tests

## Specifications

Drive of stirring points	Motorised
Stirring quantity	<b>100 ... 1000 ml</b>
Adjustable speed	0 ... 1600 U/min
Stirring surface	Ø 88 mm
Material of the platform	plastic
Power supply	Power adapter 5 V / 2 A DC 90 ... 230 V AC
Weight	ca. 500 g
Dimensions	125 x 110 x 70 mm

## More information

Manual



More product info



Similar products



Subject to change